

Technical sheet GM500 Greenline racks

Configuration	Species	Overall Dimension
cages (6 w x 9 h) Single sided	Mouse	1314 x 500 x 1742 mm
cages (6 w x 10 h) Single sided	Mouse	1314 x 500 x 1742 mm
cages (7 w x 9 h) Single sided	Mouse	1528 x 500 x 1742 mm
cages (7 w x 10 h) Single sided	Mouse	1528 x 500 x 1742 mm
cages (8 w x 10 h) Single sided	Mouse	1742 x 500 x 1742 mm
cages (8 w x 10 h) Single sided	Mouse	1742 x 500 x 1742 mm
cages (6 w x 9 h) Double sided	Mouse	1314 x 500 x 1742 mm
cages (6 w x 10 h) Double sided	Mouse	1314 x 500 x 1742 mm
cages (7 w x 9 h) Double sided	Mouse	1528 x 500 x 1742 mm
cages (7 w x 10 h) Double sided	Mouse	1528 x 500 x 1742 mm
cages (8 w x 10 h) Double sided	Mouse	1742 x 500 x 1742 mm
cages (8 w x 10 h) Double sided	Mouse	1742 x 500 x 1742 mm

€ 0,00

* with removable