

Tecniplast rabbit cage racks - Type P - 2 levels - > 5 kg



Price on demand

The P-Type cage provides a generous floor area of 5456 cm² (5.87 ft²) and has overall dimensions (W x D x H) of 870 x 740 x 131 mm (34.25 x 29.13 x 5.16 in), with an internal height of 604 mm (23.76 in). The P-Type rack holds 2 cages in a 1-wide x 2-high configuration. Its overall dimensions are 937 (970^{*}) x 845 (895^{**}) x 1718 mm, equivalent to 36.89 (38.19^{*}) x 33.27 (35.25^{**}) x 67.64 in. The cage is made of Noryl plastic. The cage body and waste tray are autoclavable and highly stackable. The waste tray is symmetric (no back , no front). Included: Perforated trays (bottom of the cage) Solid trays (to recover the excrements) Platforms (resting plate) Automatic drinking nipples Pellet feeder Hay hopper * with interlinking device ** with feeder and water bottle