

Tecniplast rabbit cage racks - Type X - 3 levels - < 5 kg



Price on demand

The X-Type cage offers a floor area of 4264 cm² (4.59 ft²) and has overall dimensions (W x D x H) of 713 x 716 x 115 mm (28.06 x 28.19 x 4.53 in), with an internal height of 476 mm (18.74 in). The X-Type rack holds 3 cages in a 1-wide x 3-high configuration. Its overall dimensions are 784 (817^{*}) x 820 (869^{**}) x 1830 mm, equivalent to 30.87 (32.17^{*}) x 32.28 (34.20^{**}) x 72.05 in. The cage is made of Noryl plastic. The cage body is autoclavable and highly stackable. The waste tray is also autoclavable and highly stackable, with a symmetric design without a fixed back or front. Included: - Perforated trays (bottom of the cage) - Solid trays (to recover the excrements) - Platforms (resting plate) - Automatic drinking nipples - Pellet feeder Option: - Hay hopper (84 pieces available) * with interlinking device ** with feeder and water bottle

Specificaties